

Latent Print Development Kits



Small Particle Reagent Kit

Use this kit to locate, develop and lift latent prints with our small particle reagent. Simply spray this special solution on a latent print and watch the print form a gray deposit that can be lifted with our lifting tape. It can be used on all nonporous surfaces, whether wet or dry! Also can be used on glass or glossy surfaces, waxed and plastic-coated paper surfaces, exterior automobile surfaces and windows.

Small Particle Reagent Kit comes with:

- 1-Plastic Carrying Case with Molded Inserts
- 2-No. 2555R1 500-ml bottles of Small Particle Reagent Dark
- 1-No. 2555W1 500-ml bottle of Small Particle Reagent Light
- 1-No. 2555UV1 500-ml bottle of Small Particle Reagent Fluorescent
- 2-No. 2556 Dispenser Spray Heads
- 2-No. 25551 Plastic 500-ml bottles for Rinse Water
- 6-No. 6705 Disposable Vinyl Gloves
- 1-No. 9696 5" Pair of Scissors
- 1-No. 7000 16-ounce aerosol can of Magic Skin Cleaner
- 1-No. 4753 2" x 360" roll of Transparent Lifting Tape
- 1-No. 4755 4" x 360" roll of Transparent Lifting Tape
- 1-No. 2100 Latent Backing Book
- 1-No. 2557 Storage Box (for towels)
- Instruction Booklet

CATALOG NO.	DESCRIPTION	PRICE
2500	Small Particle Reagent Kit	
2555R	500-ml Reagent Refill	
2560	250-ml Reagent Refill*	

*Comes with a spray head.

Adhesive Side Powder Kit

This kit contains all the components necessary to develop latent prints on the adhesive surface of tapes, labels and similar surfaces. Adhesive-side powder yields excellent results on duct tape, plastic tapes, (clear, frosted and opaque), paper labels and tapes (except those with water-activated adhesives), vinyl surgical tape (excellent fidelity but low contrast). ACE's No. ADH100 powder is ideal for vinyl or plastic black electrical tape. Use it with the ultraviolet lamp supplied. The No. ATR1 adhesive tape release agent allows tangled tape to be released without damage to latent prints which may be present. Complete instructions are included.



Adhesive Side Powder Kit comes with:

- 1-No. ADH100 Adhesive Side Powder, 50 g
- 1-No. ADH100UV Fluorescent Adhesive Side Powder, 25 g
- 1-No. ADH101 EZFLO Super Concentrate, 76 oz. (198 ml)
- 1-No. ATR1 Adhesive Tape Release Agent, 1 oz. (29 ml)
- 1-No. 9500 Longwave Ultraviolet Light, DC
- 4-AA-Cell Batteries, Alkaline
- 1-LDPE Bottle with Cap, 6 oz. (198 ml)
- 2-No. 1506 Regular Powder Brushes
- 2-No. 8607 Plastic Tweezers
- 1-No. 8823 Rinse Basin/Kit Case
- 1-No. 8825 Mixing Bowl with Lid
- 2-No. 6705 Disposable Latex Glove Pairs
- 1-No. 86051 Measuring Spoon, 1 tsp. (5 ml)
- 1-No. 86052 Measuring Spoon, 0.15 ml
- 1-No. ADH1001 Instructions for No. ASP150 Use
- 1-No. ATR11 Instructions for No. TRA20 Use

CATALOG NO.	DESCRIPTION	PRICE
ADH1000	Adhesive Side Powder Kit	

